



Mouse Anti-Feline pancreatic lipase Monoclonal Antibody, Clone 15G10 (CABT- JX53018)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	PNLIP
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Feline
Clone	15G10
Purification	Protein A purified
Conjugate	Unconjugated
Applications	<p>LFIA, CLIA, ELISA (Det)</p> <p>We recommend the following as antibody pair (Capture - Detection): CABT-JX53017-CABT-JX53019</p> <p>Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.</p>
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS with 0.1% Sodium Azide

Preservative	0.1% Sodium Azide
Storage	Store at 2-8°C, do not need to freeze.
Ship	Wet ice