



# Mouse Anti-Canine SAA Monoclonal Antibody, Clone 23F10 (CABT-JX53024)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                    |   |
|--------------------|---|
| Target             | SAA   |
| Isotype            | IgG   |
| Source/Host        | Mouse   |
| Species Reactivity | Canine  |
| Clone              | 23F10   |
| Purification       | Protein A purified  |
| Conjugate          | Unconjugated  |
| Applications       | LIA, LFIA, ELISA (Det)<br>We recommend the following as antibody pair (Capture - Detection):<br>CABT-JX53023-CABT-JX53024<br>Each laboratory should determine an optimum working titer for use in its particular application.<br>Other applications have not been tested but use in such assays should not necessarily be excluded. |
| Format             | Purified, Liquid  |
| Concentration      | Lot specific  |
| Size               | 1 mg  |
| Buffer             | PBS with 0.1% Sodium Azide  |
| Preservative       | 0.1% Sodium Azide   |

|                |  |
|----------------|--|
| <b>Storage</b> | Store at 2-8°C, do not need to freeze. |
|----------------|--|

---

|             |         |
|-------------|---------|
| <b>Ship</b> | Wet ice |
|-------------|---------|

---