



Mouse Anti-Feline AGP Monoclonal Antibody, Clone 9F21 (CABT-JX53025)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	AGP
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Feline
Clone	9F21
Purification	Protein A purified
Conjugate	Unconjugated
Applications	LIA, LFIA, ELISA (Cap) We recommend the following as antibody pair (Capture - Detection): CABT-JX53025-CABT-JX53026 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS with 0.1% Sodium Azide
Preservative	0.1% Sodium Azide

Storage	Store at 2-8°C, do not need to freeze.
----------------	--

Ship	Wet ice
-------------	---------
