



Mouse Anti-Human Calprotectin Monoclonal Antibody, Clone IN420 (CABT-JXM06C)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Native Human Calprotectin
Immunogen	Purified Human calprotectin
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	IN420
Purification	Protein A purified
Conjugate	Unconjugated
Applications	<p>ELISA, LFIA</p> <p>We recommend the following for sandwich ELISA (Capture - Detection):</p> <p>CABT-JXM06C - CABT-JXM06D</p> <p>Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.</p>
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	10 mM Phosphate Buffered Saline, pH 7.0, 0.09% Sodium Azide

Preservative	0.09% Sodium Azide
Storage	Short Term (\leq 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Ship	Wet ice