



Mouse Anti-Human Calprotectin Monoclonal Antibody, Clone IN431 (CABT-JXM06D)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Native Human Calprotectin
Immunogen	Purified Human calprotectin
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	IN431
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, LFIA We recommend the following for sandwich ELISA (Capture - Detection): CABT-JXM06C - CABT-JXM06D Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	10 mM Phosphate Buffered Saline, pH 7.0, 0.09% Sodium Azide

Preservative	0.09% Sodium Azide
Storage	Stored at +2-8°C for 36 months.
Ship	Wet ice