



Mouse Anti-Human Pancreatic Elastase Monoclonal Antibody, Clone IN2145 (CABT-JXM07C)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Human Pancreatic Elastase No reactivity with pig elastase; pig, bovine, human hemaglobin; pig, bovine, human transferrin; human calprotectin; bovine, human lactoferrin
Isotype	IgA
Source/Host	Mouse
Species Reactivity	Human
Clone	IN2145
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA (Cap) We recommend the following as antibody pair (Capture - Detection): CABT-JXM07C - CABT-JXM07D Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Liquid
Concentration	Lot specific
Size	1 mg

Buffer	10 mM Phosphate Buffered Saline, pH 7.2, 0.09% Sodium Azide
Preservative	0.09% Sodium Azide
Storage	Stored at +2-8°C
Ship	Wet ice