



# Mouse Anti-Human Haptoglobin Monoclonal Antibody, Clone IN817 (CABT-JXM08C)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	Human Haptoglobin
<b>Immunogen</b>	Purified Human Haptoglobin
<b>Isotype</b>	IgG1
<b>Source/Host</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Clone</b>	IN817
<b>Purification</b>	Protein A purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	<p>LFIA, ELISA (Cap)</p> <p>We recommend the following for sandwich ELISA (Capture - Detection): CABT-JXM08C-CABT-JXM08D</p> <p>Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.</p>
<b>Format</b>	Liquid
<b>Concentration</b>	Lot specific
<b>Size</b>	1 mg
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.0, 0.09% Sodium Azide

<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Ship</b>	Wet ice