



Mouse Anti-Human Haptoglobin Monoclonal Antibody, Clone IN818 (CABT-JXM08D)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Human Haptoglobin
Immunogen	Purified Human Haptoglobin
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	IN818
Purification	Protein A purified
Conjugate	Unconjugated
Applications	LFIA, ELISA (Det) We recommend the following for sandwich ELISA (Capture - Detection): CABT-JXM08C-CABT-JXM08D Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	10 mM Phosphate Buffered Saline, pH 7.0, 0.09% Sodium Azide

Preservative	0.09% Sodium Azide
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Ship	Wet ice