



# Rabbit Anti-Human IgA ( $\alpha 1$ & $\alpha 2$ ) Monoclonal Antibody, Clone SN239 (CABT-JXM09C)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	Reacts to human IgA, including both IgA1 and IgA2. No cross reactivity with human IgG, IgM, IgD, or IgE.
<b>Immunogen</b>	Human IgA
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Clone</b>	SN239
<b>Purification</b>	Protein A purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA (Cap), ICC, IHC We recommend the following as antibody pair (Capture - Detection): CABT-JXM09C - CABT-JXM09D Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
<b>Format</b>	Liquid
<b>Concentration</b>	1.0 mg/mL
<b>Size</b>	100 $\mu$ g
<b>Buffer</b>	50% Glycerol/PBS with 1% BSA and 0.09% sodium azide

<b>Preservative</b>	0.09% sodium azide
<b>Storage</b>	Stable for 1 Year at -20.0°C from date of receipt
<b>Ship</b>	Wet ice