



Mouse Anti-Human GPI Monoclonal Antibody, Clone 25G5 (CABT-JXM14C)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Human GPI
Isotype	IgG1
Source/Host	Mouse
Species Reactivity	Human
Clone	25G5
Purification	Protein A purified
Conjugate	Unconjugated
Applications	CLIA (Cap), ELISA (Cap) We recommend the following as antibody pair (Capture - Detection): CABT-JXM14C - CABT-JXM14D Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Preservative	None

Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Ship	Wet ice
