



Mouse anti Human Inhibin B monoclonal antibody, clone 6 (CABT-L006M)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Human Inhibin B
Immunogen	Recombinant full length protein.
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Clone	6
Purification	Protein A purified
Conjugate	Unconjugated
Applications	CG, ELISA, CLIA
Molecular Weight	150 kDa
Format	Liquid
Concentration	Lot specific
Size	0.5 mg, 1 mg
Buffer	0.01 M PBS, pH7.4
Preservative	None
Storage	Store this product at 2-8°C. Once dissolved, aliquot and store the product at -20°C. Refrain

from multiple freeze/thaw cycles.

Ship

Wet ice

BACKGROUND

Introduction

Inhibin-B is glycoproteins of the transforming growth factor- β superfamily, is produced almost exclusively by the Sertoli cells and has been proposed as direct markers of their function and indirect markers of spermatogenesis.

Keywords

Inhibin B;INHBB;INHBB;inhibin, beta B;inhibin beta B chain;Inhibin, beta-2;activin beta-B chain;activin AB beta polypeptide
