



Humanized Anti-Oxacillinase-48 type carbapenemases monoclonal antibody, clone OXA2 (CABT-L0110Y)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|----------------------|--|
| Specificity | Oxacillinase-48 type carbapenemases |
| Isotype | IgG |
| Source/Host | Humanized |
| Clone | OXA2 |
| Purification | Purity > 90% |
| Conjugate | Unconjugated |
| Applications | Suitable for use in ELISA, CLIA, and LFIA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. |
| Format | Liquid, Purified |
| Concentration | Lot specific |
| Size | 1 ea |
| Buffer | PBS or Tris |
| Storage | Store at -20°C for long-term storage. Store at 2-12°C for up to one month. Avoid freeze/thaw cycles. |
| Ship | Wet ice |