



Humanized anti Human ZnT8 monoclonal antibody, clone T2 (CABT-L0181Y)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Isotype	IgG
Source/Host	Human
Species Reactivity	Human
Clone	T2
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, CLIA Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Liquid
Concentration	Lot specific
Size	1 ea
Buffer	PBS or Tris
Preservative	None
Storage	Store at -20 °C. Stable for 12 months from date of receipt.
Ship	Wet ice