



Rabbit Anti-Human RGS5(aa 21-71) polyclonal antibody (CABT-L0198Y)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Human RGS5
Immunogen	A synthesized peptide derived from human RGS5(Accession O15539), corresponding to amino acid residues I21-N71.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purified
Conjugate	Unconjugated
Applications	Suitable for WB, IHC, IF/ICC Recommended dilution: WB 1:500-1:1000, IHC 1:50-1:200, IF/ICC 1:100-1:500 Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Liquid. Purified
Concentration	Lot specific
Size	100 ul, 200 ul
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol
Preservative	0.02% sodium azide

Storage Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Ship Wet ice

BACKGROUND

Keywords RGS5; Regulator of G-protein signaling 5; RGS5
