



Mouse Anti-Human antithrombin III monoclonal antibody, clone G5 (CABT-L0214Y)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Clone	G5
Purification	Affinity purified
Conjugate	Unconjugated
Applications	Suitable for ELISA, CLIA We recommend the following for sandwich ELISA (Capture - Detection): CABT-L0214Y - CABT-L0215Y Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Liquid
Concentration	Lot specific
Size	500 ug
Buffer	PBS or Tirs with preservative
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Ship	Wet ice