



Mouse Anti-HPV L1 Monoclonal antibody, clone 6419 (CABT-L038M)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Reactive with the L1 of types 6, 11, 16, 18, 31, 33, 45, 52, & 58.
Immunogen	Human Papillomavirus L1 capsid
Isotype	IgG
Source/Host	Mouse
Species Reactivity	HPV
Clone	6419
Purification	Affinity purification. Purity:>95% Determined by SDS-PAGE
Conjugate	Unconjugated
Applications	Suitable for use in ELISA, WB and IHC Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Buffer	0.01 M phosphate buffered saline, pH 7.2 and contains 0.1% sodium azide No stabilizing proteins.
Preservative	0.1% Sodium Azide

Storage Store at -20°C.

Ship Wet ice

BACKGROUND

Keywords HPV; Human papillomavirus;
