



Rabbit Anti-Human DDX19A monoclonal Antibody, clone 341343E4 (CABT-L0412Y)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Human DDX19A
Immunogen	Human DDX19A fusion protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Clone	341343E4.
Purification	Affinity purified
Conjugate	Unconjugated
Applications	Suitable for use in Indirect ELISA, ELISA(Cap), and Cytometric Bead Array We recommend the following for sandwich ELISA (Capture - Detection): CABT-L0412Y - CABT-L0413Y Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg
Buffer	PBS

Storage	Upon receipt, prepare single use aliquots and store at -20 °C until use. Avoid multiple freeze thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	ATP-dependent RNA helicase involved in mRNA export from the nucleus.
Keywords	DDX19 DDX19L, DDX19 like protein, DDX19A, DDX19L, DEAD box protein 19A, DKFZp686C21137, FLJ11126

GENE INFORMATION

Entrez Gene ID	55308
UniProt ID	Q9NUU7