



Human LECT2 ELISA Matched Antibody Pair Set (CABT-L0598Y)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	<p>Capture antibody: CABT-C0598Y Unconjugated Recombinant Rabbit Anti Human LECT2 monoclonal antibody (IgG) supplied in PBS-only buffer, 1 vial, 100ug (~ 1 mg/ml) Detection antibody: CABT-D0598Y Unconjugated Recombinant Rabbit Anti Human LECT2 monoclonal antibody (IgG) supplied in PBS-only buffer, 1 vial, 100ug (~ 1 mg/ml)</p>
Specificity	Human LECT2
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity Purification
Conjugate	Unconjugated
Applications	<p>Sandwich ELISA</p> <p>Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.</p>
Format	Purified, Liquid
Concentration	Lot specific
Size	1 set (Reagents are sufficient for 10 x 96 well plates)
Buffer	PBS

Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Ship	Wet ice

BACKGROUND

Introduction	Has a neutrophil chemotactic activity. Also a positive regulator of chondrocyte proliferation.
Keywords	chm II; chm2; Chondromodulin II; hLECT2; LECT 2; LECT-2; LECT2; LECT2_HUMAN; Leukocyte cell derived chemotaxin 2; Leukocyte cell-derived chemotaxin-2; MGC126628;

GENE INFORMATION

Entrez Gene ID	3950
UniProt ID	O14960