



Rabbit Anti Human LECT2 monoclonal Antibody, Clone 8T6 (CABT-L0599Y)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Recombinant Rabbit anti Human LECT2 Monoclonal Antibody
Specificity	Human LECT2
Immunogen	Recombinant protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Clone	8T6
Purification	Affinity Purification
Conjugate	Unconjugated
Applications	sELISA, ELISA(Cap) We recommend the following for sandwich ELISA (Capture - Detection): CABT-L0599Y - CABT-L0600Y Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Size	100 µg

Buffer	0.01M PBS, pH 7.2
Storage	For short-term storage, store at 2-8°C. For long-term storage, aliquot and store at -20°C. Avoid multiple freeze-thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	Has a neutrophil chemotactic activity. Also a positive regulator of chondrocyte proliferation (PubMed:9524238). Research Field: Signal Transduction
Keywords	LECT2; hLECT2

GENE INFORMATION

Entrez Gene ID	3950
UniProt ID	O14960