



Rabbit anti Cat Fel d 4 polyclonal antibody (CABT-L075)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	The antiserum contains IgG antibodies to Fel d 4 as well as IgG antibodies to other cat allergens.
Target	Fel d 4
Immunogen	Fel d 4
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Cat
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	100 µl
Buffer	1%BSA 50%glycerol PBS, pH 7.4.
Preservative	None
Storage	Store at 4°C.

BACKGROUND

Introduction Cat allergy is unique among allergy to mammals in that the allergen Fel d 4 is a uteroglobin-like

protein and not a lipocalin. The biochemical spectrum of the cat allergens is thus uncertain, particularly with regard to the role that a cat lipocalin protein may play in sensitization to cats in allergic individuals.

Keywords	Allergen Fel d 4; uteroglobin-like protein;
-----------------	---------------------------------------------

GENE INFORMATION

UniProt ID	Q5VFH6
-------------------	------------------------
