



Rabbit anti Guinea Pig Guinea Pig Urinary Protein polyclonal antibody (CABT-L076)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Binds to guinea pig urinary proteins. Reacts with guinea pig epithelial and hair extracts. Does not recognize rat or mouse urinary proteins, does not react with rabbit or hamster epithelial and hair extracts.
Target	Urinary protein
Immunogen	Guinea pig urinary protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Guinea pig
Purification	Affinity chromatography using recombinant Protein G.
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Concentration	1-2 mg/ml
Size	100 µl
Buffer	1-2mg/ml in phosphate buffered saline, pH 7.4.
Preservative	None
Storage	Store at 4°C.

BACKGROUND

Introduction	Major urinary proteins (Mups), also known as α 2u-globulins, are a subfamily of proteins found in abundance in the urine and other secretions of many animals.
Keywords	Major urinary protein; Mups; α 2u-globulins