



# Rabbit anti Arachis hypogaea Ara h 2 polyclonal antibody (CABT-L083)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                    |                                                  |
|--------------------|--------------------------------------------------|
| Specificity        | The antiserum contains IgG antibodies to Ara h 2 |
| Target             | Ara h 2                                          |
| Immunogen          | Natural Ara h 2                                  |
| Isotype            | IgG                                              |
| Source/Host        | Rabbit                                           |
| Species Reactivity | Arachis hypogaea                                 |
| Conjugate          | Unconjugated                                     |
| Applications       | ELISA                                            |
| Format             | Liquid                                           |
| Size               | 100 µl                                           |
| Buffer             | PBS, pH 7.4, containing 1%BSA/50% glycerol.      |
| Preservative       | None                                             |
| Storage            | Store at 4°C.                                    |

## BACKGROUND

|              |                                                                    |
|--------------|--------------------------------------------------------------------|
| Introduction | Ara h 2 is a seed storage protein from Arachis hypogaea (peanuts). |
|--------------|--------------------------------------------------------------------|

**Keywords**

Arachis hypogaea; Conglutin; 2S protein 1; Seed storage protein SSP1; Seed storage protein SSP2; IgE-binding protein; Protease inhibitor; Seed storage; Allergen Ara h 2; Peanuts

---

## GENE INFORMATION

**UniProt ID**

[Q6PSU2](#)

---