



Rabbit Anti Glycyphagus domesticus Gly d 2 polyclonal antibody (CABT-L086)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Binds to species specific epitopes present on Glycyphagus domesticus allergen, Gly d 2.
Target	Gly d 2
Immunogen	Gly d 2
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Glycyphagus domesticus
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	100 µl
Buffer	PBS, pH 7.2.
Preservative	None
Storage	Store at 4°C.

BACKGROUND

Introduction	Of 17 G.domesticus positive patients, 16 showed IgE binding to rGly d 2 on immunoblot.
--------------	--

Keywords

Mite group 2 allergen Gly d 2.02; Mite; Glycyphagus domesticus

GENE INFORMATION

UniProt ID

[Q9U5P7](#)
