



Rabbit Anti Lepidoglyphus destructor Lep d 2 polyclonal antibody (CABT-L087)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Binds to species specific epitopes present on Lepidoglyphus destructor allergen, Lep d 2.
Target	Lep d 2
Immunogen	Lep d 2
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Lepidoglyphus destructor
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	100 µl
Buffer	PBS, pH 7.2.
Preservative	None
Storage	Store at 4°C.

BACKGROUND

Keywords	Lep d 2; Lep d 1
----------	------------------

GENE INFORMATION

UniProt ID

[P80384](#)
