



Sheep Anti-Human Hyaluronic acid polyclonal antibody (CABT-L1000)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Sheep anti hyaluronic acid antibody recognizes hyaluronic acid, a mucopolysaccharide which is an integral part of the gel-like substance of animal connective tissue that serves as a lubricating and protective agent.
Target	Hyaluronic acid
Immunogen	Hyaluronic acid from human umbilical cord
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Broad
Purification	Purified IgG prepared by affinity chromatography on Protein G
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Concentration	5 mg/ml
Size	1 ml
Buffer	Phosphate buffered saline
Preservative	0.09% Sodium Azide
Storage	Store at +4°C or at -20°C if preferred. Storage in frost-free freezers is not recommended. This

product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.

BACKGROUND

Introduction

Hyaluronic acid is a mucopolysaccharide which is an integral part of the gel-like substance of animal connective tissue that serves as a lubricating and protective agent.

Keywords

Hyaluronan;hyaluronate;HA;Hyaluronic Acid
