



Mouse Anti FLAG monoclonal antibody, clone N3 (CABT-L1031)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	This affinity purified monoclonal antibody binds to fusion proteins containing a FLAG peptide sequence. The antibody recognizes the FLAG peptide sequence at the N-terminus, Met-N-terminus, C-terminus, and internal sites of the fusion protein.
Specificity	Binding site: N-Asp-Tyr-Lys-Asp-Asp-Asp-Lys-C
Target	FLAG
Immunogen	FLAG; peptide sequence DYKDDDDK
Isotype	IgG1
Source/Host	Mouse
Clone	N3
Purification	Affinity purified
Conjugate	Unconjugated
Applications	WB, IP, IHC, IF, ICC
Preparation	Preparation Note: Immediately before use, dilute in Tris buffered saline (TBS), pH 8.0, with 3% nonfat Milk
Format	Liquid
Concentration	1 mg/ml
Size	200 µg, 1 mg, 5 mg

Buffer	Contains 50% glycerol in phosphate buffered saline.
Preservative	None
Storage	-20° C, Avoid Freeze/Thaw Cycles
Ship	Wet ice

BACKGROUND

Keywords Flag (DYKDDDDK), FITC-labeled;DDDDK epitope tag;DYKDDDDK epitope tag;DYKDDDDK tag;ECS epitope tag;ECS tag;Enterokinase Cleavage Site epitope tag;Enterokinase Cleavage Site tag;flag;FLAG;Flag