



Rabbit Anti-Human ZNF652 polyclonal antibody (CABT-L1054)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	ZNF652
Immunogen	Synthetic peptide directed towards the N-terminal of human ZNF652
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Bovine, Dog, Horse, Goat, Human, Mouse, Pig, Rabbit, Rat, Yeast
Purification	Affinity Chromatography
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Concentration	0.5 mg/ml
Size	100 µl
Buffer	PBS with 2% sucrose
Preservative	0.09% Sodium Azide
Storage	-20° C, Avoid Freeze/Thaw Cycles

BACKGROUND

Introduction	ZNF652 is a new candidate transcription factor.
Keywords	ZNF652;zinc finger protein 652;KIAA0924;DKFZp781E2122;KIAA0924;OTTHUMP00000220047;Zinc finger protein 652;ZN652_HUMAN;ZNF652

GENE INFORMATION

Entrez Gene ID	22834
UniProt ID	Q9Y2D9