



# Rabbit Anti Mouse GPR68 polyclonal antibody (CABT-L1058)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Target</b>	GPR68
<b>Immunogen</b>	This antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 200-225 amino acids from the Central region of mouse Gpr68.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Mouse
<b>Purification</b>	Protein A purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Format</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze-thaw cycles.

## BACKGROUND

<b>Keywords</b>	GPR68;G protein-coupled receptor 68;ovarian cancer G-protein coupled receptor 1;OGR1;G protein coupled receptor 68;G-protein coupled receptor 68;GPR 12A;GPR 68;GPR12A;Gpr68
-----------------	--