



Mouse Anti-Human lambda light chain monoclonal antibody, clone 4E23 (CABT-L1089)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------------|---|
| Specificity | Specific for independent constantly expressed epitope of lambda chain. Reacts with monomer and dimer forms of free lambda chains but not with lambda chain-bearing immunoglobulins. |
| Immunogen | Mixture of light chains of human immunoglobulin. |
| Isotype | IgG2a |
| Source/Host | Mouse |
| Species Reactivity | Human |
| Clone | 4E23 |
| Purification | Protein A purified |
| Conjugate | Unconjugated |
| Applications | WB, IHC, IHC-Fr, IHC-P |
| Format | Liquid |
| Size | 200 µg |
| Buffer | PBS (pH 7.4) |
| Preservative | 0.09% Sodium Azide |
| Storage | 4° C. Do not freeze. |

BACKGROUND

Keywords

IGLC;immunoglobulin lambda constant group
