



# Rabbit Anti-Human CXCL9 polyclonal antibody (CABT-L1096)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Target	CXCL9
Immunogen	KLH conjugated synthetic peptide between 96-122 amino acids from the C-terminal region of human CXCL9
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity chromatography
Conjugate	Unconjugated
Applications	IHC, WB
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	-20° C, Avoid Freeze/Thaw Cycles

## BACKGROUND

Introduction	The function of CXCL9 has not been specifically defined; however, it is thought to be involved in T cell trafficking. This protein has been localized to 4q21 with INP10, which is also a member of the chemokine family of cytokines.
--------------	--

**Keywords**

CXCL9;chemokine (C-X-C motif) ligand 9;CMK;MIG;Humig;SCYB9;crg-10;C-X-C motif  
chemokine 9;small-inducible cytokine B9;gamma-interferon-induced monokine;monokine  
induced by gamma interferon;monokine induced by interferon-gamma

---

## GENE INFORMATION

**Entrez Gene ID**[4283](#)**UniProt ID**[Q07325](#)