



Rabbit Anti-CHD1 polyclonal antibody (CABT-L1134)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	CHD1
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	ELISA, WB, IP, IHC, IF
Molecular Weight	200 kDa
Format	Liquid
Size	150 µl
Buffer	PBS with 50% glycerol pH 7.3.
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquoting is unnecessary for -20°C storage.

BACKGROUND

Introduction	The CHD family of proteins is characterized by the presence of chromo (chromatin organization modifier) domains and SNF2-related helicase/ATPase domains. CHD genes alter gene expression possibly by modification of chromatin structure thus altering access of the transcriptional apparatus to its chromosomal DNA template.
Keywords	CHD1;chromodomain helicase DNA binding protein 1;chromodomain-helicase-DNA-binding protein 1;CHD-1;ATP-dependent helicase CHD1;DKFZp686E2337

GENE INFORMATION

Official Symbol	CHD1
Synonyms	ATP dependent helicase CHD1, CHD 1, CHD1, DKFZp686E2337
Entrez Gene ID	1105
UniProt ID	B3KT33