



RHA™ Mouse Anti-Iprodione monoclonal antibody, clone IPD (CABT-L1148)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Iprodione [KLH]
Isotype	IgG
Source/Host	Mouse
Clone	IPD
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS, pH 7.4
Preservative	None
Storage	-20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Iprodione is used on crops affected by Botrytis bunch rot, Brown rot, Sclerotinia and other fungal diseases in plants. It is currently applied in a variety of crops: fruit, vegetables, ornamental trees
--------------	---

and scrubs and on lawns. It is a contact fungicide that inhibits the germination of fungal spores and it blocks the growth of the fungal mycelium.

Keywords

Iprodione
