



RHA™ Mouse Anti-Spinosyn monoclonal antibody, clone SPS (CABT-L1149)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Spinosyn [KLH]
Isotype	IgG
Source/Host	Mouse
Clone	SPS
Purification	Purified from mouse ascites.
Conjugate	Unconjugated
Applications	ELISA, LFIA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS, pH 7.4
Preservative	None
Storage	-20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords	Spinosyn
----------	----------