



Rabbit Anti-Mouse AKT Pan monoclonal antibody, clone K.425.5 (CABT-L1185)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	AKT1
Immunogen	Synthetic peptide corresponding to residues in the carboxy-terminal sequence of mouse Akt
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Fruit fly, Human, Mouse, Non-human primate, Pig, Rat
Clone	K.425.5
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	FC, ICC, IHC-P, IP, IF, WB
Format	Liquid
Buffer	0.01M HEPES, pH 7.5, with 0.15M NaCl, 100µg/ml BSA, 50% glycerol
Preservative	See individual product datasheet
Storage	-20°C

BACKGROUND

Keywords AKT1;v-akt murine thymoma viral oncogene homolog 1;RAC-alpha serine/threonine-protein kinase;AKT;PKB;PRKBA;RAC;PKB alpha;RAC-PK-alpha;proto-oncogene c-Akt;protein kinase

GENE INFORMATION

Entrez Gene ID [11651](#)

UniProt ID [P31750](#)
