



Rabbit Anti-Mouse BAD (phospho Ser112) monoclonal antibody, clone H.556.0 (CABT-L1207)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	BAD
Immunogen	Synthetic phosphopeptide corresponding to residues surrounding pSer112 of mouse Bad
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Non-human primate, Rat
Clone	H.556.0
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	FC, IHC-P, WB
Format	Liquid
Buffer	0.01M HEPES, pH 7.5, with 0.15M NaCl, 100µg/ml BSA, 50% glycerol
Preservative	See individual product datasheet
Storage	-20°C

BACKGROUND

Keywords

BAD;BCL2-associated agonist of cell death;BBC2;BCL2L8;bcl2-associated agonist of cell death;bcl2-L-8;BCL2-binding protein;bcl-2-like protein 8;BCL2-binding component 6;bcl-2-binding component 6

GENE INFORMATION

Entrez Gene ID[12015](#)**UniProt ID**[Q61337](#)