



Rabbit Anti-Human Torc1 monoclonal antibody, clone U.303.3 (CABT-L1252)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	CRTC1
Immunogen	Synthetic peptide corresponding to residues surrounding Val463 of human TORC1 (CRTC1) protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Clone	U.303.3
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	0.01M HEPES, pH 7.5, with 0.15M NaCl, 100µg/ml BSA, 50% glycerol
Preservative	See individual product datasheet
Storage	-20°C

BACKGROUND

Keywords CRTC1;CREB regulated transcription coactivator 1;MECT1;TORC1;TORC-1;WAMTP1;CREB-

regulated transcription coactivator 1;mucoepidermoid carcinoma translocated protein
1;transducer of regulated cAMP response element-binding protein (CREB) 1

GENE INFORMATION

Entrez Gene ID [23373](#)

UniProt ID [Q6UUV9](#)
