



Rabbit Anti-Human Torc3 monoclonal antibody, clone S.296.5 (CABT-L1253)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	CRTC3
Immunogen	Synthetic peptide corresponding to the sequence of human TORC3 (CRTC3) protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Non-human primate, Rat
Clone	S.296.5
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	0.01M HEPES, pH 7.5, with 0.15M NaCl, 100µg/ml BSA, 50% glycerol
Preservative	See individual product datasheet
Storage	-20°C

BACKGROUND

Keywords	CRTC3;CREB regulated transcription coactivator 3;CREB-regulated transcription coactivator 3;FLJ21868;TORC-3;transducer of CREB protein 3;transducer of regulated CREB protein
-----------------	---

GENE INFORMATION

Entrez Gene ID [64784](#)

UniProt ID [Q6UUV7](#)
