



Rabbit Anti-Human CSK monoclonal antibody, clone D.789.5 (CABT-L1254)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody is not cross-reactive with other proteins.
Target	CSK
Immunogen	Synthetic peptide corresponding to residues surrounding amino acid residue Val399 of human Csk
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Non-human primate, Pig, Rat
Clone	D.789.5
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	0.01M HEPES, pH 7.5, with 0.15M NaCl, 100µg/ml BSA, 50% glycerol
Preservative	See individual product datasheet
Storage	-20°C

BACKGROUND

Keywords

CSK;c-src tyrosine kinase;tyrosine-protein kinase CSK;C-Src kinase;protein-tyrosine kinase
CYL

GENE INFORMATION

Entrez Gene ID[1445](#)**UniProt ID**[P41240](#)