



Rabbit Anti-Human HSP40 monoclonal antibody, clone H.880.6 (CABT-L1269)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody is not cross-reactive with other HSPs.
Target	DNAJB1
Immunogen	Synthetic peptide corresponding to Glu223 of human HSP40/Hdj1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Non-human primate, Rat
Clone	H.880.6
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	IHC-P, WB
Format	Liquid
Buffer	0.01M HEPES, pH 7.5, with 0.15M NaCl, 100µg/ml BSA, 50% glycerol
Preservative	See individual product datasheet
Storage	-20°C

BACKGROUND

Keywords

DNAJB1;DnaJ (Hsp40) homolog, subfamily B, member

1;Hdj1;Sis1;HSPF1;Hsp40;RSPH16B;dnaJ homolog subfamily B member 1;hDj-1;human DnaJ protein 1

GENE INFORMATION

Entrez Gene ID[3337](#)**UniProt ID**[P25685](#)
