



Rabbit Anti-Human EEA1 monoclonal antibody, clone G.54.2 (CABT-L1273)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	EEA1
Immunogen	Synthetic peptide corresponding to residues surrounding pSer70 of human EEA1 protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Clone	G.54.2
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	ICC, IP, IF, WB
Format	Liquid
Buffer	0.01M HEPES, pH 7.5, with 0.15M NaCl, 100µg/ml BSA, 50% glycerol
Preservative	See individual product datasheet
Storage	-20°C

BACKGROUND

Keywords	EEA1;early endosome antigen 1;MST105;ZFYVE2;MSTP105;early endosome antigen 1, 162kD;endosome-associated protein p162;early endosome-associated protein;zinc finger
----------	--

GENE INFORMATION

Entrez Gene ID [8411](#)

UniProt ID [Q15075](#)
