



Sheep Anti-Norsertraline Polyclonal antibody (CABT-L128S)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Norsertraline antibody was raised in sheep using Norsertraline-BTG as the immunogen
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	Ig fraction supplied in liquid form in 20 mM phosphate buffer, 150 mM sodium chloride and 0.09% sodium azide, pH 7.2
Preservative	0.09% Sodium Azide
Storage	Maintain at -20°C for up to 12 months. Avoid repeated freeze-thaw cycles. Store product undiluted.
Ship	Wet ice

BACKGROUND

Introduction	Norsertraline is an antidepressant of the selective serotonin reuptake inhibitor (SSRI) class.
Keywords	Zoloft; Sertraline; Norsertraline
