



Rabbit Anti-Human FOXD3 monoclonal antibody, clone C.92.9 (CABT-L1308)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	FOXD3
Immunogen	Synthetic peptide corresponding to residues at the C-terminus of human FoxD3
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Clone	C.92.9
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	0.01M HEPES, pH 7.5, with 0.15M NaCl, 100µg/ml BSA, 50% glycerol
Preservative	See individual product datasheet
Storage	-20°C

BACKGROUND

Introduction	This gene belongs to the forkhead family of transcription factors which is characterized by a distinct forkhead domain. Mutations in this gene cause autoimmune susceptibility 1.
--------------	---

Keywords FOXD3;forkhead box D3;forkhead box protein D3;Genesis;HNF2;AIS1;forkhead box D3;Forkhead box protein D3;foxd3;FOXD3_HUMAN;Genesis;HNF2;HNF3/FH transcription factor genesis;OTTHUMP00000010713;HNF3/FH transcription factor genesis;AIS1

GENE INFORMATION

Entrez Gene ID [27022](#)

UniProt ID [Q9UJU5](#)
