



Rabbit Anti-Human GCN5 monoclonal antibody, clone T.886.5 (CABT-L1317)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This antibody is not cross-reactive with the related PCAF protein.
Target	GCN5
Immunogen	Synthetic peptide corresponding to the amino terminus of human GCN5L2
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Bovine, Horse, Human, Mouse, Non-human primate, Rat
Clone	T.886.5
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	ICC, IP, IF, WB
Format	Liquid
Buffer	0.01M HEPES, pH 7.5, with 0.15M NaCl, 100µg/ml BSA, 50% glycerol
Preservative	See individual product datasheet
Storage	-20°C

BACKGROUND

Introduction	KAT2A, or GCN5, is a histone acetyltransferase.
Keywords	GCN5;Gcn5p;AAS104;YGR252W;ADA4;SWI9;NP_011768.1;EC 2.3.1.48;Eukaryota;Fungi;Ascomycota;Saccharomycotina (true yeasts);Taphrinomycotina;Schizosaccharomycetes (fission yeasts);Basidiomycota;Agaricomycotina;Tremellomycetes;Pucciniomycotina;Microbotryomycetes

GENE INFORMATION

Entrez Gene ID	2648
-----------------------	----------------------

UniProt ID	Q92830
-------------------	------------------------