



Mouse Anti-ssDNA/dsDNA monoclonal Antibody, clone N4E9 (CABT-L139M)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Binds unspecifically to double- and single-stranded DNA oligomers. This antibody binds to DNA oligomers both single- and double-stranded, irrespective of the sequence, showing DNA-hydrolytic activity.
Immunogen	3'-biotin oligodeoxynucleotide ss(dN)40 with N=CCATGAGTGATAACACTGCGGCCAACTTACTTCTGACAAC
Isotype	IgG1, κ
Source/Host	Mouse
Species Reactivity	N/A
Clone	N4E9
Purification	Affinity purification
Conjugate	Unconjugated
Applications	Suitable for use in hyrdolysis, SPR, and ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Purified, Liquid
Concentration	Lot specific
Buffer	PBS with 0.02% Proclin 300.
Storage	Store at -20°C.

Ship

Wet ice

BACKGROUND

Keywords

Single Stranded DNA; ssDNA; DNA (Single Stranded); DNA
