



Rabbit Anti-Human JMJD1B monoclonal antibody, clone T.404.0 (CABT-L1438)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	KDM3B
Immunogen	Synthetic peptide corresponding to the sequence of the human JMJD1B protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Non-human primate
Clone	T.404.0
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	ICC, IP, IF, WB
Format	Liquid
Buffer	0.01M HEPES, pH 7.5, with 0.15M NaCl, 100µg/ml BSA, 50% glycerol
Preservative	See individual product datasheet
Storage	-20°C

BACKGROUND

Keywords KDM3B;lysine (K)-specific demethylase 3B;5qNCA;NET22;C5orf7;JMJD1B;lysine-specific demethylase 3B;nuclear protein 5qNCA;jumonji domain containing 1B;jumonji domain-

GENE INFORMATION

Entrez Gene ID [51780](#)

UniProt ID [Q7LBC6](#)
