



Rabbit Anti-LY6E polyclonal antibody (CABT-L144R)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	LY6E
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-101 of human LY6E (NP_002337.1).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Purification	Affinity purification
Conjugate	Unconjugated
Applications	WB, IF
Sequence Similarities	LMCFSCLNQKSNLYCLKPTICSDQDNYCVTVSASAGIGNLVTFGHSLSKTCSPACPIPEG VNVGVASMGISCCQSFLCNFS
Molecular Weight	Observed MW: 14kDa; Calculated MW:13kDa
Cellular Localization	Cell membrane, GPI-anchor, Lipid-anchor
Positive Control	LO2, 293T, SKOV3, MCF7, U-87MG, Mouse brain, Mouse kidney, Mouse liver
Format	Liquid
Concentration	Lot specific

Size	50 µl
Buffer	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Preservative	0.02% sodium azide
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Ship	Wet ice

GENE INFORMATION

Entrez Gene ID [4061](#)