



Rabbit Anti-Human MEK1/MEK2 monoclonal antibody, clone F.002.8 (CABT-L1453)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	MAP2K1; MAP2K2
Immunogen	Synthetic peptide corresponding to the sequence of human MEK1/2
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Dog, Human, Mouse, Non-human primate, Rat
Clone	F.002.8
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	0.01M HEPES, pH 7.5, with 0.15M NaCl, 100µg/ml BSA, 50% glycerol
Preservative	See individual product datasheet
Storage	-20°C

BACKGROUND

Keywords MAP2K1;MAP2K2;MAPKK1;MAPKK2;MEK1;MEK2;MEKK1;MK2;MKK1;MKK2;PRKMK1;PRKMK2