



Rabbit Anti-Human MESDC2 monoclonal antibody, clone T.640.0 (CABT-L1471)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	MESDC2
Immunogen	Recombinant human MESDC2 protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Non-human primate, Rat
Clone	T.640.0
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	IP, WB
Format	Liquid
Buffer	0.01M HEPES, pH 7.5, with 0.15M NaCl, 100µg/ml BSA, 50% glycerol
Preservative	See individual product datasheet
Storage	-20°C

BACKGROUND

Keywords	MESDC2;mesoderm development candidate 2;LDLR chaperone MESD;BOCA;KIAA0081;MESD;mesoderm development protein;renal carcinoma antigen NY-
----------	--

GENE INFORMATION

Entrez Gene ID [23184](#)

UniProt ID [Q14696](#)
