



Rabbit Anti-Human MNK1 monoclonal antibody, clone F.722.5 (CABT-L1476)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	MKNK1
Immunogen	Synthetic peptide corresponding to the sequence of human Mnk1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Clone	F.722.5
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	0.01M HEPES, pH 7.5, with 0.15M NaCl, 100µg/ml BSA, 50% glycerol
Preservative	See individual product datasheet
Storage	-20°C

BACKGROUND

Keywords MKNK1;MAP kinase interacting serine/threonine kinase 1;MNK1;MAP kinase-interacting serine/threonine-protein kinase 1;MAPK signal-integrating kinase 1;MAP kinase signal-

GENE INFORMATION

Entrez Gene ID	8569
----------------	----------------------

UniProt ID	Q9BUB5
------------	------------------------
