



Rabbit Anti-Human GBL monoclonal antibody, clone F.270.3 (CABT-L1479)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	MLST8
Immunogen	Synthetic peptide corresponding to residues surrounding Gln210 of human GBL
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Non-human primate, Rat
Clone	F.270.3
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	IP, WB
Format	Liquid
Buffer	0.01M HEPES, pH 7.5, with 0.15M NaCl, 100µg/ml BSA, 50% glycerol
Preservative	See individual product datasheet
Storage	-20°C

BACKGROUND

Keywords	MLST8;MTOR associated protein, LST8 homolog (S. cerevisiae);GBL;LST8;POP3;WAT1;GbetaL;target of rapamycin complex subunit
-----------------	---

GENE INFORMATION

Entrez Gene ID [64223](#)

UniProt ID [Q9BVC4](#)
